



TGEX-FH Expression Vector

INSTRUCTION MANUAL

TGEX-FH Transient Mammalian Expression Vector

Catalog #: MX004

Version: A1.2 – September 2020

Table of Contents

Limited Use License for the TGEX Vector Series	4
Description	5
Introduction	5
Content, Shipping & Storage	5
Limited Product Warranty	5
TGEX Vector Series	6
Vector Map	6
Cloning Site	7
Feature Table	7
Restriction Site Summary	8
Experimental Procedures	10
General Molecular Biology Techniques	10
Plasmid Maintenance	10
Cloning into TGEX- FH	10
Sequencing of Inserts	12
Fab Fragment Expression	12
Appendix	13
MSDS Information	13
Quality Control	13
Technical Support	13
References	13

Limited Use License for the TGEX Vector Series

As a condition of sale of this product to you, and prior to using this product, you must agree to the terms and conditions of this license. Antibody Design Labs grants to the buyer with the sale of any of its **TGEX™** vectors (the “Product”) a non-exclusive, non-transferable and limited license to use the Product in research only conducted by the buyer. Such license specifically excludes the right to sell or otherwise transfer the Product, its components or derivatives thereof to third parties. No modifications to the Product may be made without express written permission from Antibody Design Labs. The buyer is not granted a license to use the Product for human or animal therapeutic, diagnostic, or prophylactic purposes.

Antibody Design Labs does not warrant that the use or sale of the Product, the use thereof in combination with other products, or the use of the Product in the operation of any process will not infringe the claims of any United States or other patent(s).

If the buyer is not willing to accept the limitations of this license, without modification, buyer may refuse this license by returning the Product unopened and unused. By keeping or using the Product, buyer agrees to be bound by the terms of this license.

Entities wishing to use the Product for commercial purposes are required to obtain a license from Antibody Design Labs. Commercial purposes may include, but are not limited to: use of the Product in manufacturing, use of the Product to provide a service, use of the Product for therapeutic or diagnostic purposes, or resale of the Product, whether or not such Product is resold for use in research. For information on purchasing a commercial license to the Product, please contact a licensing representative by phone at (858) 480-6213 or by e-mail at licensing@abdesignlabs.com.

All trademarks are the property of their respective owners.

Description

Introduction

The **TGEX™** vector series is designed for the rapid expression of antibody molecules by transient gene expression in mammalian cell suspension culture. The vector series features a CMV promoter, the adenovirus tripartite leader sequence (TPL) (Logan 1984, Mariati 2010), a variable antibody domain leader sequence with its intron and convenient cloning sites to insert antibody variable regions. The constant regions are derived from the human IgG1 and kappa sequences. There is no selection marker and the overall small size of the vectors is optimized for transient transfection; expensive antibiotics to prepare large quantities of plasmid for transient transfections are not required. Using widely available cell lines and large scale transfection technologies (see experimental procedure) yields of antibody between 10 mg/L and 100 mg/L in serum-free conditions are routinely achieved in the laboratory in just a few days.

The **TGEX™-FH** vector is designed for the expression of a heavy chain variable region with the CH1 constant region of the human IgG1 heavy chain. Expression of Fab fragment molecules is achieved by co-transfection with a light chain variable region cloned into the vector **TGEX™-LC**.

Content, Shipping & Storage

Content

VECTOR	COMPOSITION	AMOUNT
TGEX-FH	20 µl at 0.5 µg/µl of DNA vector in DNA Conservation Buffer (Tris-HCL 5 mM, EDTA 0.1 mM, pH 8.5)	10 µg

Shipping & Storage

TGEX-FH vector is shipped on wet ice. Upon receipt, store the vector at -20°C.

Limited Product Warranty

This warranty limits our liability to the replacement of this product. No other warranties of any kind express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Antibody Design Labs. Antibody Design Labs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

For research use only; not intended for any animal or human therapeutic or diagnostic use.

TGEX Vector Series

Characteristics of the TGEX vector series

VECTOR	CATALOG	USE	DESCRIPTION
TGEX™-HC	MX001	Cloning of VH domain	For the expression of human IgG1 heavy chain
TGEX™-LC	MX002	Cloning of VL domain	For the expression of human Kappa light chain
TGEX™-FC	MX003	Cloning of VH domain	For the expression of human IgG1 Fc fusion.
TGEX™-FH	MX004	Cloning of VH domain	For the expression of human IgG1 Fab fragment
TGEX™-HChis	MX005	Cloning of VH domain	For the expression of human IgG1 heavy chain with a HIS tag
TGEX™-SCblue	MX006	scFv cloning vector	For the transfer of scFv from any PADL phagemid vector and expression as an scFv-Fc fusion


Combination of vectors to desired antibody format

COMBINATION	FORMAT	PURIFICATION ¹
TGEX™-HC + TGEX™-LC	Full length human or chimeric human IgG1/K	Protein A or G
TGEX™-HChis + TGEX™-LC	Full length human or chimeric human IgG1/K	Protein A , G or IMAC
TGEX™-FH + TGEX™-LC	Human or chimeric IgG1/K Fab fragment	Protein L, G or IMAC
TGEX™-FC	Fc fusion	Protein A or G
TGEX™-SCblue	scFv-Fc fusion	Protein A or G

1. Purification by protein A, G and L may require testing.


Vector Map

The figure below illustrates the main features of **TGEX-FH** expression vector. The full vector sequence is available online for download in varied formats on the product web page; the total length of the vector is 3668 bp.



Cloning Site

Following is an illustration of **TGEX-FH** cloning site from the EcoRI site and onward. The VH domain is inserted between the BssHII site and the BsmBI site; after excision of the intron sequence, the IGHV1-2 sequence encodes the leader peptide MDWTWRILFLVAAATGAHS. SA/SD donor and acceptor sites.



Feature Table

The features of **TGEX-FH** transient expression vector are highlighted in the following table.

FEATURE	LOCATION	DESCRIPTION
Promoter	5-585	CMV promoter.
TPL	612-1000	Adenovirus tripartite leader sequence (Logan 1984, Mariati 2010)
IGHV1-2 leader	1025-1192	Human IGHV1-2 leader sequence with intron. The mature RNA encodes the 18 a.a.-long signal sequence MDWTWRILFLVAAATGAHS; cleavage occurs on the C-terminal side of the terminal serine.
Human G1 CH1	1209-1556	Sequence encoding the human IgG1 heavy chain sequence comprising the CH1, CH2 and CH3 domains with the hinge region. The sequence is intronless; the C-terminal lysine residue is encoded.
HIS tag	1533-1550	6x HIS tag located before a C-terminal glycine residue
BGpA	2225 -1678	Rabbit beta-globin polyadenylation signal sequence.
pMB1 origin	1752- 2371	pBR322 origin for replication in <i>E. coli</i> with a temperature-sensitive high copy-number phenotype (Lin-Chao 1992).
TEM1 beta-lactamase	2526- 3386	Ampicillin resistance for selection in <i>E. coli</i> .

Restriction Site Summary

Enzyme	Site	Nb	Position	Strand	Isoschizomers
AgeI	A [^] CCGGT	1	1313		AsiGI BshTI CspAI PinAI
AlwNI	CAGNNN [^] CTG	1	2117		CaiI PstNI
ApaI	GGGCC [^] C	1	1220		Bsp120I PspOMI
ArsI	(8/13)GACNNNNNNNTTYG (11/6)	1	851		
BalI	TGG [^] CCA	1	3596		MlsI MluNI MscI Msp20I
BamHI	G [^] GATCC	1	1524		
BbsI	GAAGAC (2/6)	1	1230	-	BpiI BpuAI BstV2I
BbvCI	CCTCAGC (-5/-2)	1	1400		
BcgI	(10/12)CGANNNNNNNTGC (12/10)	1	3115		
BdaI	(10/12)TGANNNNNNNTCA (12/10)	1	1788		
BlpI	GC [^] TNAGC	1	1552		Bpu1102I Bsp1720I CelIII
BsePI	G [^] CGCGC	1	1182		BssHII PauI PteI
Bsp1407I	T [^] GTACA	1	3630		BsrGI BstAUI
BstEII	G [^] GTNACC	1	1412		BstPI Eco91I EcoO65I PspEI
Bsu36I	CC [^] TNAGG	1	1384		Eco81I AxyI Bse21I
DrdI	GACNNNN [^] NNGTC	1	1808		AasI DseDI
EagI	C [^] GGCCG	1	1558		BseX3I BstZI EclXI Eco52I
Eam1105I	GACNNN [^] NNGTC	1	2594		AhdI AspEI BmeRI DriI
EcoNI	CCTNN [^] NNNAGG	1	1289		BstENI XagI
EcoRI	G [^] AATTC	1	1019		
Esp3I	CGTCTC (1/5)	1	1202		BsmBI
FspI	TGC [^] GCA	1	2819		Acc16I AviII NsbI
HindIII	A [^] AGCTT	1	1701		
NarI	GG [^] CGCC	1	1338		DinI EgeI EheI KasI Mly113I SfoI SspDI
NmeAIII	GCCGAG (21/19)	1	2722		
NotI	GC [^] GGCCGC	1	1557		CciNI
PasI	CC [^] CWGGG	1	1279		
PmeI	GTTT [^] AAAC	1	1000		MssI
PvuI	CGAT [^] CG	1	2966		Ple19I BpvUI MvrI
PvuII	CAG [^] CTG	1	642		
SacII	CCGC [^] GG	1	740		Sfr303I KspI SgrBI Cfr42I SstII
SnaBI	TAC [^] GTA	1	357		BstSNI Eco105I
SpeI	A [^] CTAGT	1	18		AhlI BcuI
Tth111I	GACN [^] NNGTC	1	1319		AspI PflFI PsyI
XhoI	C [^] TCGAG	1	966		StrI TliI Sfr274I PaeR7I SlaI
XmnI	GAANN [^] NNTTC	1	3194		Asp700I MroXI PdmI
AclI	AA [^] CGTT	2	2824		Psp1406I
		2	3197		
BciVI	GTATCC (6/5)	2	1915	-	BfuI BsuI
		2	3442		
Bpu10I	CCTNAGC (-5/-2)	2	930		
		2	1400		
BseRI	GAGGAG (10/8)	2	715		
		2	1247		
BseYI	CCCAGC (-5/-1)	2	1473		GsaI
		2	2010		
BspHI	T [^] CATGA	2	2426		CciI PagI RcaI
		2	3434		
BsrBI	CCGCTC (-3/-3)	2	1555	-	AccBSI MbiI
		2	3438		
BssSI	CACGAG (-5/-1)	2	1879	-	BauI Bst2BI
		2	3263		
BstXI	CCANNNNN [^] NTGG	2	1426		

		2	1541	
BtgZI	GCGATG (10/14)	2	368	-
		2	623	
BtsI	GCAGTG (2/0)	2	2992	
		2	3020	
CspCI	(11/13) CAANNNNNGTGG (12/10)	2	407	
		2	1508	
FalI	(8/13) AAGNNNNNCTT (13/8)	2	784	
		2	1296	
NcoI	C^CATGG	2	379	Bsp19I
		2	1049	
NdeI	CA^TATG	2	252	FauNDI
		2	1657	
OliI	CACNN^NNGTG	2	1160	AleI
		2	1481	
PflMI	CCANNNN^NTGG	2	641	AccB7I BasI Van91I
		2	1049	
PpiI	(7/12) GAACNNNNNCTC (13/8)	2	2397	
		2	3240	
SacI	GAGCT^C	2	583	Ecl136II EcoICRI Eco53kI
				Psp124BI SstI
		2	1034	
VspI	AT^TAAT	2	25	AseI PshBI
		2	2770	

Absent Sites:

AarI, AbsI, AflII, AjuI, AlfI, AloI, AscI, AsuII, AvrII, BaeI, BarI, BclI, BplI, BsaBI, BsaXI, BsgI, BsiWI, BsmI, BspEI, BspMI, BstAPI, BstZ17I, BtrI, ClaI, DraIII, Eco31I, Eco47III, EcoRV, FseI, FspAI, HpaI, KflI, KpnI, MauBI, MfeI, MluI, MreI, NaeI, NheI, NruI, NsiI, PacI, PciI, PfoI, PmaCI, PshAI, PsiI, PspXI, PsrI, PstI, RsrII, SalI, SapI, SexAI, SfiI, SgfI, SgrAI, SgrDI, SmaI, SphI, SrfI, Sse8387I, StuI, SwaI, TstI, XbaI, XcmI.

Experimental Procedures

General Molecular Biology Techniques

Molecular biology should be conducted under the supervision of a qualified instructor trained to standard safety practice in a molecular biology laboratory environment. Standard molecular biology procedures can be found in a general molecular biology handbook such as Sambrook (1989).

Plasmid Maintenance

Propagation and maintenance of **TGEX-FH** is obtained on any *recA1*, *endA1* *E. coli* strain using LB or 2xYT medium supplemented with ampicillin (100 µg/ml) as a selection marker and incubated at 37°C with agitation. **TGEX-FH** is a derivative of pBR322 with a high copy number origin of replication and usually gives high yields of plasmid DNA with most standard laboratory strains such as XL1-blue or DH5α. The high copy number phenotype is temperature sensitive and requires incubation at 37°C (Lin-Chao 1992). Some DNA stabilizing strains are known to produce smaller amounts of plasmid DNA. In case of issues, we recommend using XL10-Gold® from Agilent Technologies, Inc., on which **TGEX-FH** plasmid DNA can be isolated in large quantities.

Cloning into TGEX-FH

In Silico Clone Design

A complete IGHV1-2 leader sequence MDWTWRILFLVAAATGAHS is necessary for secretion of the antibody in the culture supernatant and proper removal of the leader peptide by host proteases. In the following schema, after cutting by BssHII and BsmBI, the end of the leader encoding sequence 5'-CGCGCACTCC will be removed. This short sequence must be included back in the final clone for proper protein maturation and secretion. The BsmBI site will be eliminated during the cloning.

```

                                     BssHII                               human G1 CH1
                                     |                               |-----|
                                     AlaHisSer                     AlaSerThrLysGlyProSer
1161  ACTTCTGTGTTCTCTCCACAGGCGCGCACTCC [VH DOMAIN] GCATCCACCAAGGGCCCATCG

```

Vector Digestion

BssHII and BsmBI restriction enzymes are fully active at 50°C and 55°C, respectively. The double digestion of the vector can be either conducted at 50°C with both enzymes together or at 55°C with BsmBI initially alone and together with BssHII after cooling the reaction to 50°C. Esp3I is an isoshizomer of BsmBI with an optimal working temperature of 37°C. Please, consult the documentation of your restriction enzyme provider for optimal double digestion conditions.

Primer Design for Restriction Cloning with the BssHII Site

Oligo1 is an example of primer designed to amplify a VH domain sequence and clone it into the BssHII site. A minimum of 2 nucleotides is recommended to cut BssHII site close to the end (source New England's BioLabs); these 2 N are followed by the BssHII site GCGCGC and the end of the leader sequence ACTCC encoding the end of the leader peptide; the resulting NNGCGCGCACTCC extension is added 5' to the VH domain primer.

Oligo1 5' -NNGCGCGCACTCC-VH-start

Primer Design for Restriction Cloning with the BsmBI Site

Oligo2 is an example of primer designed to amplify the VH domain from the end of the J region and includes a BsmBI site compatible with **TGEX-FH** cloning site. BsmBI is a type IIS restriction enzyme that cuts outside of its recognition site. BsmBI will cut immediately before the start of the human IgG1 constant region sequence, exactly 1 base after the end of the site and 5 bases further on the opposite strand, thus freeing a 5' 4-bases overhang GCAT on the sense strand. In Oligo2, a BsmBI site situated symmetrically to the vector will generate a complementary overhang. After digestion and ligation, the two BsmBI sites, the one in the vector and the one in Oligo2, will be removed. We added two nucleotides on the 5' end of the primer although a minimum of one nucleotide is recommended to cut BsmBI site close to the end (source New England BioLabs).

Oligo2 5' -NNcgtctcNATGC-JH-end

Removing BsmBI Sites in J Regions

As shown in Table 1, many heavy chain J regions contain a BsmBI site near their 3' end. These sites can be eliminated by modifying the codon usage in the sequence encoding the TVS amino acid sequence. We suggest using GTT instead of GTC for the codon encoding the valine residue; the resulting mismatch will create minimal perturbation and simply adding one extra nucleotide to the J primer is sufficient to ensure strong amplification of the VH domain.

J Region	3' Sequence
Human IGHJ1	ctggtcac cgctctc ctca g
Human IGHJ2	ctggtcactgtctctctca g
Human IGHJ3	atggtcac cgctctc ttca g
Human IGHJ4	ctggtcac cgctctc ctca g
Human IGHJ5	ctggtcac cgctctc ctca g
Human IGHJ6	acggtcac cgctctc ctca g
Mouse mus. IGHJ1*01	acggtcac cgctctc ctca g
Mouse mus. IGHJ1*02	acggtcaccggtttctca g
Mouse mus. IGHJ1*03	acggtcac cgctctc ctca g
Mouse mus. IGHJ2*01	actctcacagtctctctca g
Mouse mus. IGHJ2*02	tctctcacagtctctctca g
Mouse mus. IGHJ2*03	agtctcacagtctctctca g
Mouse mus. IGHJ3*01	ctggtcactgtctctctgca g
Mouse mus. IGHJ4*01	tcagtcac cgctctc ctca g
Rabbit IGHJ1	ctggtcaccatctcttca g
Rabbit IGHJ2	ctggtcac cgctctc ctca g
Rabbit IGHJ3	ctggtcac cgctctc ctca g
Rabbit IGHJ4	ctggtcac cgctctc ctca g
Rabbit IGHJ5*01	ctggtcac cgctctc ttca n
Rabbit IGHJ5*02	ctggtcactgtctcttca g
Rabbit IGHJ6	ctcgtcac cgctctc ttca g
Translation	T V S S/A

Table 1. Occurrence of BsmBI sites on the 3' end of human, mouse musculus and rabbit functional J regions (source IMGT⁵).

Alternative to BsmBI Sites

All restriction enzymes that generate 4-bases long 5' overhangs can be used in place of BsmBI; this is the case for example of BsaI (GGTCTC(1/5)), another type IIS restriction enzyme. This option could come in handy when the VH domain contains another BsmBI preventing cloning.

Sequencing of Inserts

The following primers give a strong PCR amplification of the TGEX vector series inserts between the EcoRI site and the NotI site. The primer *tgex-S3* can be used to sequence the VH domain in full.

tgex-S3 5'- AGGCGTCTAACCAGTCACAGTC

tgex-R2 5'- CAAAAAATCCAACACACTATTGC

Fab Fragment Expression

Cell lines

Cell lines adapted for culture in suspension and serum-free conditions are recommended. HEK293 and CHO cells are often used for antibody expression by transient transfection; you can either adapt your own cell line or get it from a supplier (e.g. Life Technology). HEK293 cells are particularly well suited for expression using **TGEX™** vector series.

Transient Transfection

Many transfection reagents especially designed for transient transfection are commercially available from multiple providers (e.g. Life Technologies, Mirus Bio LLC). We recommend testing the transfection conditions with a reporter plasmid first to determine the percentage of cells effectively transfected and optimal transfection conditions; fluorescent reporters are often used with that purpose. Similarly, any condition known to boost expression should be carefully tested in your system before being scaled up. We did observe a strong increase in expression in 293 cells upon exposure to sodium valproate (Backliwal 2008).

LIGHT CHAIN TO HEAVY CHAIN RATIO

We recommend starting with a 1:1 light chain to heavy chain ratio during transfection. We observed many antibodies with a better expression at a 2:1 light chain to heavy chain ratio although each antibody requires fine tuning for optimal expression.

IMAC Purification

Follow the recommendation of your IMAC resin provider for purification using the HIS tag. The HIS tag located at the C-terminal of the heavy chain can sustain wash in the presence of up to 10 mM imidazole at pH 7.4; leakiness may appear at 20 mM and above.

Appendix

MSDS Information

MSDSs (Material Safety Data Sheets) are available on **Antibody Design Labs** website at the corresponding product page.

Quality Control

Specifications and quality control are detailed on the online product page. **Antibody Design Labs** certifies that the product will perform according to these specifications.

Technical Support

Visit **Antibody Design Labs** website at www.abdesignlabs.com for technical resources, including manuals, vector maps and sequences, application notes, FAQs, etc.

For more information or technical assistance, call, write, fax, or email us at:

Antibody Design Labs

4901 Morena Blvd, Suite 203
San Diego, CA 92117

Email: support@abdesignlabs.com

Phone: 1-877-223-3104 (Toll Free)

Fax: 1-858-272-6007 (24 hour)

(Monday – Friday 9:00 AM – 5:00 PM PST)

References

1. LOGAN, J., & SHENK, T. (1984). ADENOVIRUS TRIPARTITE LEADER SEQUENCE ENHANCES TRANSLATION OF MRNAS LATE AFTER INFECTION. *PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES OF THE UNITED STATES OF AMERICA*, 81(12), 3655–9.
2. MARIATI, HO, S. C. L., YAP, M. G. S., & YANG, Y. (2010). EVALUATING POST-TRANSCRIPTIONAL REGULATORY ELEMENTS FOR ENHANCING TRANSIENT GENE EXPRESSION LEVELS IN CHO K1 AND HEK293 CELLS. *PROTEIN EXPRESSION AND PURIFICATION*, 69(1), 9–15.
3. LIN-CHAO S, CHEN WT, WONG TT (1992). HIGH COPY NUMBER OF THE PUC PLASMID RESULTS FROM A ROM/ROP-SUPPRESSIBLE POINT MUTATION IN RNA II. *MOL MICROBIOL*. 22:3385-93.
4. SAMBROOK, J., FRITSCH, E.F., AND MANIATIS, T. (1989). IN MOLECULAR CLONING: A LABORATORY MANUAL. COLD SPRING HARBOR LABORATORY PRESS, NY, VOL. 1, 2, 3.
5. IMGT®, THE INTERNATIONAL IMMUNOGENETICS INFORMATION SYSTEM® [HTTP://WWW.IMGT.ORG](http://www.imgt.org) (FOUNDER AND DIRECTOR: MARIE-PAULE LEFRANC, MONTPELLIER, FRANCE).
6. BACKLIWAL G, HILDINGER M, KUETTEL I, DELEGRANGE F, HACKER DL, WURM FM. (2008). VALPROIC ACID: A VIABLE ALTERNATIVE TO SODIUM BUTYRATE FOR ENHANCING PROTEIN EXPRESSION IN MAMMALIAN CELL CULTURES. *BIOOTECHNOL BIOENG*, 101(1):182-9.

This product is subject to Antibody Design Labs Terms & Conditions of Sales available online at <http://www.abdesignlabs.com/terms/>.

© 2015 Antibody Design Labs. All rights reserved.
