



TGEX-LC-hL Expression Vector

INSTRUCTION MANUAL

TGEX-LC-hL Transient Mammalian Expression Vector

Catalog #: MX013

Version: A1.2 – August 2020

Table of Contents

Limited Use License for the TGEX Vector Series	4
Description	5
Introduction	5
Content, Shipping & Storage	5
Limited Product Warranty	5
TGEX Vector Series	6
Vector Map	7
Cloning Site	7
Feature Table	7
Restriction Site Summary	8
Experimental Procedures	10
General Molecular Biology Techniques	10
Plasmid Maintenance	10
Cloning into TGEX-LC-hL	10
Sequencing of Inserts	11
Antibody Expression	11
Appendix	12
MSDS Information	12
Quality Control	12
Technical Support	12
References	12

Limited Use License for the TGEX Vector Series

As a condition of sale of this product to you, and prior to using this product, you must agree to the terms and conditions of this license. Antibody Design Labs grants to the buyer with the sale of any of its **TGEX™** vectors (the “Product”) a non-exclusive, non-transferable and limited license to use the Product in research only conducted by the buyer. Such license specifically excludes the right to sell or otherwise transfer the Product, its components or derivatives thereof to third parties. No modifications to the Product may be made without express written permission from Antibody Design Labs. The buyer is not granted a license to use the Product for human or animal therapeutic, diagnostic, or prophylactic purposes.

Antibody Design Labs does not warrant that the use or sale of the Product, the use thereof in combination with other products, or the use of the Product in the operation of any process will not infringe the claims of any United States or other patent(s).

If the buyer is not willing to accept the limitations of this license, without modification, buyer may refuse this license by returning the Product unopened and unused. By keeping or using the Product, buyer agrees to be bound by the terms of this license.

Entities wishing to use the Product for commercial purposes are required to obtain a license from Antibody Design Labs. Commercial purposes may include, but are not limited to: use of the Product in manufacturing, use of the Product to provide a service, use of the Product for therapeutic or diagnostic purposes, or resale of the Product, whether or not such Product is resold for use in research. For information on purchasing a commercial license to the Product, please contact a licensing representative by phone at (858) 480-6213 or by e-mail at licensing@abdesignlabs.com.

All trademarks are the property of their respective owners.

Description

Introduction

The **TGEX™** vector series is designed for the rapid expression of antibody molecules by transient gene expression in mammalian cell suspension culture. The vector series features a CMV promoter, the adenovirus tripartite leader sequence (TPL) (Logan 1984, Mariati 2010), a variable antibody domain leader sequence with its intron and convenient cloning sites to insert antibody variable regions. The constant regions are derived from the human IgG1 and kappa sequences. There is no selection marker and the overall small size of the vectors is optimized for transient transfection; expensive antibiotics to prepare large quantities of plasmid for transient transfections are not required. Using widely available cell lines and large scale transfection technologies (see experimental procedure) yields of antibody between 10 mg/L and 100 mg/L in serum-free conditions are routinely achieved in the laboratory in just a few days.

The **TGEX™-LC-hL** vector is designed for the expression of a light chain variable region with the constant region of the human lambda light chain. Expression of full length antibody molecules is achieved by co-transfection with a heavy chain variable region cloned into one vector **TGEX™-HC** series.

Content, Shipping & Storage

Content

VECTOR	COMPOSITION	AMOUNT
TGEX-LC-hL	20 µl at 0.5 µg/µl of DNA vector in DNA Conservation Buffer (Tris-HCL 5 mM, EDTA 0.1 mM, pH 8.5)	10 µg

Shipping & Storage

TGEX-LC vector is shipped on wet ice. Upon receipt, store the vector at -20°C.

Limited Product Warranty

This warranty limits our liability to the replacement of this product. No other warranties of any kind express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Antibody Design Labs. Antibody Design Labs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

For research use only; not intended for any animal or human therapeutic or diagnostic use.

TGEX Vector Series

Characteristics of the TGEX vector series

VECTOR	CATALOG	USE	DESCRIPTION
TGEX™-HC	MX001	Cloning of VH domain	For the expression of human IgG1 heavy chain
TGEX™-LC	MX002	Cloning of VL domain	For the expression of human Kappa light chain
TGEX™-FC	MX003	Cloning of VH domain	For the expression of human IgG1 Fc fusion.
TGEX™-FH	MX004	Cloning of VH domain	For the expression of human IgG1 Fab fragment
TGEX™-HChis	MX005	Cloning of VH domain	For the expression of human IgG1 heavy chain with a HIS tag
TGEX™-SCblue	MX006	scFv cloning vector	For the transfer of scFv from any PADL phagemid vector and expression as an scFv-Fc fusion

Combination of vectors to desired antibody format

COMBINATION	FORMAT	PURIFICATION ¹
TGEX™-HC + TGEX™-LC	Full length human or chimeric human IgG1/K	Protein A or G
TGEX™-HChis + TGEX™-LC	Full length human or chimeric human IgG1/K	Protein A , G or IMAC
TGEX™-FH + TGEX™-LC	Human or chimeric IgG1/K Fab fragment	Protein L, G or IMAC
TGEX™-FC	Fc fusion	Protein A or G
TGEX™-SCblue	scFv-Fc fusion	Protein A or G


1. Purification by protein A, G and L may require testing.

TGEX vector series for desired isotype

VECTOR	CATALOG	ISOTYPE	DESCRIPTION
TGEX™-HC-hG1	MX001	Human IgG1	Heavy chain expression plasmid for human IgG1, same as TGEX-HC
TGEX™-HC-hG2	MX014	Human IgG1	Heavy chain expression plasmid for human IgG2
TGEX™-LC-hK	MX002	Human Kappa	Light chain expression plasmid for human Kappa, same as TGEX-LC
TGEX™-LC-hL	MX013	Human Lambda	Light chain expression plasmid for human Lambda

Vector Map

The figure on the right illustrates the main features of **TGEX-LC-hL** expression vector. The full vector sequence is available online for download in varied formats on the product web page; the total length of the vector is 3787 bp.



Cloning Site

Following is an illustration of **TGEX-LC** cloning site from the EcoRI site and onward. The VL domain is inserted between the BspEI site and the BsaI site; after excision of the intron sequence, the IGKV3-20 sequence encodes the leader peptide METPAQLLFLLLLWLPDTSG. SA/SD donor and acceptor sites.

```

          tgex-S3                                     EcoRI
          |----->                                     |
971      GAAAGGCGTC TAACCAGTCA CAGTCGCAAG TTTAAACGGA TCTCTAGCGA
          |
          Start of IGKV3-20 leader sequence
          |
          1043
          MetGluTh rProAlaGln LeuLeuPheL
1021      ATTCTAGGAC CCAGAGGGAA CCATGGAAAC CCCAGCGCAG CTTCTCTTCC

          euLeuleuLe uTrpLeuPro SD
1071      TCCTGCTACT CTGGCTCCCA GGTGAGGGGA ACATGGGATG GTTTTGCATG
1121      TCAGTGAAAA CCCTCTCAAG TCCTGTTACC TGGCAACTCT GCTCAGTCAA
1171      TACAATAATT AAAGCTCAAT ATAAAGCAAT AATTCTGGCT CTTCTGGGAA
1221      GACAATGGGT TTGATTTAGA TTACATGGGT GACTTTTCTG TTTTATTTCC

          BspEI  VL start  BsaI  huCL constant region
          |      |          |          |
          SA  ThrSerGly          GlyGlnProL
1271      AATCTCAGAT ACCTCCGGAG AACGTGAGTA GACGGTCTCG GGTCAGCCCA

          ysAlaAlaPr oSerValThr LeuPheProP roSerSerGl uGluLeuGln
1321      AGGCTGCCCC CTCGGTCACT CTGTTCCCGC CCTCCTCTGA GGAGCTTCAA
  
```

Feature Table

The features of **TGEX-LC** transient expression vector are highlighted in the following table.

FEATURE	LOCATION	DESCRIPTION
Promoter	5-585	CMV promoter.
TPL	612-1000	Adenovirus tripartite leader sequence (Logan 1984, Mariati 2010)
IGVK3-20 leader	1025-1289	Human IGKV3-20 leader sequence with intron. The mature RNA encodes the 20 a.a.-long signal sequence METPAQLLFLLWLPTSG; cleavage occurs on the C-terminal side of the terminal glycine.
Human CK CDS	1311-1631	Sequence encoding the human kappa light chain constant region sequence.
BGpA	1696 - 1794	Rabbit beta-globin polyadenylation signal sequence.
pMB1 origin	1868 - 2487	pBR322 origin for replication in <i>E. coli</i> with a temperature-sensitive high copy-number phenotype (Lin-Chao 1992).
TEM1 beta-lactamase	2642 - 3502	Ampicillin resistance for selection in <i>E. coli</i> .

Restriction Site Summary

Enzyme	Site	Nb	Position	Strand	Isoschizomers
AlwNI	CAGNNN^CTG	1	2233		CaiI PstNI
BalI	TGG^CCA	1	3712		MlsI MluNI Mox20I MscI Msp20I
BcgI	(10/12) CGANNNNNNTGC (12/10)	1	3231		
BsaXI	(9/12) ACNNNNNCTCC (10/7)	1	1078		
Bsp1407I	T^GTACA	1	3746		BsrGI BstAUI
BspEI	T^CCGGA	1	1284		Aor13HI Bsp13I Kpn2I BseAI MroI AccIII
BsrBI	CCGCTC (-3/-3)	1	3554	-	AccBSI MbiI
CspCI	(11/13) CAANNNNNGTGG (12/10)	1	407		
DrdI	GACNNNN^NNGTC	1	1924		AasI DseDI
EagI	C^GGCCG	1	1674		BseX3I BstZI EclXI Eco52I
Eam1105I	GACNNN^NNGTC	1	2710		AhdI BmeRI DriI
EcoRI	G^AATTC	1	1019		
FalI	(8/13) AAGNNNNNCTT (13/8)	1	784		
FspI	TGC^GCA	1	2935		Acc16I NsbI
GsuI	CTGGAG (16/14)	1	2800	-	BpmI
HindIII	A^AGCTT	1	1817		
NmeAIII	GCCGAG (21/19)	1	2838		
NotI	GC^GGCCGC	1	1673		CciNI
OliI	CACNN^NNGTG	1	1382		AleI
PflMI	CCANNNN^NTGG	1	641		AccB7I Van91I
PmeI	GTTT^AAAC	1	1000		MssI
PstI	CTGCA^G	1	1639		BspMAI
PvuI	CGAT^CG	1	3082		Ple19I
SacI	GAGCT^C	1	583		Ecl136II EcoICRI Eco53kI Psp124BI SstI
SacII	CCGC^GG	1	740		Sfr303I KspI SgrBI Cfr42I
SapI	GCTCTTC (1/4)	1	1208		BspQI LguI PciSI
SmaI	CCC^GGG	1	1412		Cfr9I TspMI XmaI
SnaBI	TAC^GTA	1	357		BstSNI Eco105I
SpeI	A^CTAGT	1	18		AhlI BcuI
Sse8387I	CCTGCA^GG	1	1638		SbfI SdaI
XhoI	C^TCGAG	1	966		Sfr274I Paer7I SlaI
XmnI	GAANN^NNTTC	1	3310		Asp700I MroXI PdmI

AclI	AA^CGTT	2	2940	Psp1406I
		2	3313	
ApaLI	G^TGCAC	2	2136	Alw44I VneI
		2	3382	
ArsI	(8/13) GACNNNNNNTTYG (11/6)	2	851	
		2	1221	
BbsI	GAAGAC (2/6)	2	1218	BpiI BstV2I
		2	1601	
BciVI	GTATCC (6/5)	2	2031	- BfuI BsuI
		2	3558	
BglIII	A^GATCT	2	1811	
		2	3635	
BmrI	ACTGGG (5/4)	2	313	- BmuI
		2	2760	
Bpu10I	CCTNAGC (-5/-2)	2	930	
		2	1536	
BseYI	CCCAGC (-5/-1)	2	1051	GsaI PspFI
		2	2126	
BspHI	T^CATGA	2	2542	CciI PagI
		2	3550	
BssSI	CACGAG (-5/-1)	2	1995	- BauI Bst2BI
		2	3379	
BtgZI	GCGATG (10/14)	2	368	-
		2	623	
Eco31I	GGTCTC (1/5)	2	1304	Bso31I BsaI BspTNI
		2	1470	
NcoI	C^CATGG	2	379	Bsp19I
		2	1041	
NdeI	CA^TATG	2	252	FauNDI
		2	1773	
PvuII	CAG^CTG	2	642	
		2	1565	
VspI	AT^TAAT	2	25	AseI PshBI
		2	2886	

Absent Sites:

AanI, AarI, AbsI, Acc36I, Acc65I, AcvI, AdeI, AfeI, AflIII, AgeI, AjiI, AjuI, AlfI, AloI, Aor51HI, ApaI, AscI, AsiGI, AsiSI, Asp718I, AspA2I, AsuII, AsuNHI, AvrII, AxyI, BaeI, BamHI, BarI, BbrPI, BbvCI, BclI, BfrI, BfuAI, BlnI, BlpI, BmgBI, BmtI, BoxI, BplI, Bpu1102I, Bpu14I, Bsa29I, BsaBI, Bse21I, Bse8I, BseCI, BseJI, BsePI, BsgI, BshTI, BshVI, BsiWI, BsmBI, BsmI, Bsp119I, Bsp120I, Bsp1720I, Bsp68I, BspDI, BspMI, BspOI, BspT104I, BspTI, BssHII, BssNAI, Bst1107I, BstAFI, BstAPI, BstBI, BstEII, BstENI, BstPAI, BstPI, BstXI, BstZ17I, Bsu15I, Bsu36I, BsuTUI, BtrI, BtuMI, BveI, ClaI, CpoI, CsiI, CspAI, CspI, DinI, DraIII, Eco147I, Eco32I, Eco47III, Eco72I, Eco81I, Eco91I, EcoNI, EcoO65I, EcoRV, EcoT22I, EgeI, EheI, Esp3I, FbaI, FseI, FspAI, HpaI, I-CeuI, I-PpoI, I-SceI, KasI, KflI, KpnI, KroI, Ksp22I, KspAI, MabI, MauBI, MfeI, MluI, Mly113I, Mph1103I, MreI, MroNI, MspCI, MunI, Mva1269I, NaeI, NarI, NgoMIV, NheI, NruI, NsiI, NspV, PI-PspI, PI-SceI, PacI, PaeI, PalAI, PasI, PauI, PceI, PciI, PctI, PdiI, Pfl23II, PflFI, PfoI, PinAI, PluTI, PmaCI, PmlI, PscI, PshAI, PsiI, PspCI, PspEI, PspLI, PspOMI, PspXI, PsrI, PsyI, PteI, RgaI, RigI, RruI, Rsr2I, RsrII, SalI, SexAI, SfaAI, SfiI, SfoI, SfuI, SgfI, SgrAI, SgrDI, SgsI, SmiI, SphI, SrfI, SseBI, SspDI, StuI, SwaI, Tth111I, Vha464I, XagI, XbaI, XcmI, XmaJI, Zsp2I.

Experimental Procedures

General Molecular Biology Techniques

Molecular biology should be conducted under the supervision of a qualified instructor trained to standard safety practice in a molecular biology laboratory environment. Standard molecular biology procedures can be found in a general molecular biology handbook such as Sambrook (1989).

Plasmid Maintenance

Propagation and maintenance of **TGEX-LC-hL** is obtained on any *recA1*, *endA1* *E. coli* strain using LB or 2xYT medium supplemented with ampicillin (100 µg/ml) as a selection marker and incubated at 37°C with agitation. **TGEX-LC-hL** is a derivative of pBR322 with a high copy number origin of replication and usually gives high yields of plasmid DNA with most standard laboratory strains such as XL1-blue or DH5α. The high copy number phenotype is temperature sensitive and requires incubation at 37°C (Lin-Chao 1992). Some DNA stabilizing strains are known to produce smaller amounts of plasmid DNA. In case of issues, we recommend using XL10-Gold® from Agilent Technologies, Inc., on which **TGEX-LC-hL** plasmid DNA can be isolated in large quantities.

Cloning into TGEX-LC-hL

In Silico Clone Design

A complete IGKV3-20 leader sequence METPAQLLFLLWLPTSG is necessary for secretion of the antibody in the culture supernatant and proper removal of the leader peptide by host proteases. In the following schema, after cutting by BspEI and BsaI, the end of the BspEI site 5'-CCGGA will be removed. The BspEI site must be restored in the final clone for proper protein maturation and secretion. The BsaI site will be eliminated during the cloning. The VL domain will be included in-frame between the last A of the leader peptide and the first C of the human CK constant region.

		<i>BspEI</i>			<i>huCL constant region</i>

		ThrSerGly			GlyGlnProL
1271	AATCTCAGAT	ACC	TCCGGA	[VL DOMAIN]	GGTCAGCCCA

Vector Digestion

BspEI and BsaI restriction enzymes are fully active at 37°C, respectively. Please, consult the documentation of your restriction enzyme provider for buffer and optimal double digestion conditions.

Primer Design for Restriction Cloning with the BspEI Site

Oligo1 is an example of primer designed to amplify a VL domain sequence and clone it into the BspEI site. Six extra-nucleotides were added before the BspEI site to ensure proper digestion close to the end; these 6 N are followed by the BspEI site TCCGGA; the resulting NNNNNNTCCGGA extension is added 5' to the VL domain primer.

Oligo1 5'-NNNNNNTCCGGA-VL-start

Primer Design for Restriction Cloning with the BsaI Site

Oligo2 is an example of primer designed to amplify the VL domain from the end of the J region and includes a BsaI site compatible with **TGEX-LC** cloning site. BsaI is a type IIS restriction enzyme that cuts outside of its recognition site. BsaI will cut immediately before the start of the human kappa constant region sequence, exactly 1 base after the end of the site and 5 bases further on the opposite strand, thus freeing a 5' 4-bases overhang CGAA on the sense strand. In Oligo2, a BsaI site situated symmetrically to the vector will generate a complementary overhang. After digestion and ligation, the two BsaI sites, the one in the vector and the one in Oligo2, will be removed. We added two nucleotides on the 5' end of the primer although a minimum of one nucleotide is recommended to cut BsaI site close to the end (source New England's BioLabs).

Oligo2 5' -NNgggtctcNTTCG-JH-end

Alternative to BsaI Sites

All restriction enzymes that generate 4-bases long 5' overhangs can be used in place of BsaI; this is the case for example of BsmBI (CGTCTC(1/5)), another type IIS restriction enzyme. This option could come in handy when the VL domain contains another BsaI preventing cloning.

Sequencing of Inserts

The following primers give a strong PCR amplification of the **TGEX™** vector series inserts between the EcoRI site and the NotI site. The primer *tgex-S3* can be used to sequence the VH domain in full.

tgex-S3 5'- AGGCGTCTAACCAGTCACAGTC

tgex-R2 5'- CAAAAAATTCCAACACACTATTGC

Antibody Expression

Cell lines

Cell lines adapted for culture in suspension and serum-free conditions are recommended. HEK293 and CHO cells are often used for antibody expression by transient transfection; you can either adapt your own cell line or get it from a supplier (e.g. Life Technology). HEK293 cells are particularly well suited for expression using **TGEX™** vector series.

Transient Transfection

Many transfection reagents especially designed for transient transfection are commercially available from multiple providers (e.g. Life Technologies, Mirus Bio LLC). We recommend testing the transfection conditions with a reporter plasmid first to determine the percentage of cells effectively transfected and optimal transfection conditions; fluorescent reporters are often used with that purpose. Similarly, any condition known to boost expression should be carefully tested in your system before being scaled up. We did observe a strong increase in expression in 293 cells upon exposure to sodium valproate (Backliwal 2008).

LIGHT CHAIN TO HEAVY CHAIN RATIO

We recommend starting with a 1:1 light chain to heavy chain ratio during transfection. We observed many antibodies with a better expression at a 2:1 light chain to heavy chain ratio although each antibody requires fine tuning for optimal expression.

Appendix

MSDS Information

MSDSs (Material Safety Data Sheets) are available on **Antibody Design Labs** website at the corresponding product page.

Quality Control

Specifications and quality control are detailed on the online product page. **Antibody Design Labs** certifies that the product will perform according to these specifications.

Technical Support

Visit **Antibody Design Labs** website at www.abdesignlabs.com for technical resources, including manuals, vector maps and sequences, application notes, FAQs, etc.

For more information or technical assistance, call, write, fax, or email us at:

Antibody Design Labs

4901 Morena Blvd, Suite 203

San Diego, CA 92117

Email: support@abdesignlabs.com

Phone: 1-877-223-3104 (Toll Free)

Fax: 1-858-272-6007 (24 hour)

(Monday – Friday 9:00 AM – 5:00 PM PST)

References

1. LOGAN, J., & SHENK, T. (1984). ADENOVIRUS TRIPARTITE LEADER SEQUENCE ENHANCES TRANSLATION OF MRNAS LATE AFTER INFECTION. *PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES OF THE UNITED STATES OF AMERICA*, 81(12), 3655–9.
2. MARIATI, HO, S. C. L., YAP, M. G. S., & YANG, Y. (2010). EVALUATING POST-TRANSCRIPTIONAL REGULATORY ELEMENTS FOR ENHANCING TRANSIENT GENE EXPRESSION LEVELS IN CHO K1 AND HEK293 CELLS. *PROTEIN EXPRESSION AND PURIFICATION*, 69(1), 9–15.
3. LIN-CHAO S, CHEN WT, WONG TT (1992). HIGH COPY NUMBER OF THE PUC PLASMID RESULTS FROM A ROM/ROP-SUPPRESSIBLE POINT MUTATION IN RNA II. *MOL MICROBIOL*. 22:3385-93.
4. SAMBROOK, J., FRITSCH, E.F., AND MANIATIS, T. (1989). IN MOLECULAR CLONING: A LABORATORY MANUAL. COLD SPRING HARBOR LABORATORY PRESS, NY, VOL. 1, 2, 3.
5. BACKLIWAL G, HILDINGER M, KUETTEL I, DELEGRANGE F, HACKER DL, WURM FM. (2008). VALPROIC ACID: A VIABLE ALTERNATIVE TO SODIUM BUTYRATE FOR ENHANCING PROTEIN EXPRESSION IN MAMMALIAN CELL CULTURES. *BIOOTECHNOL BIOENG*, 101(1):182-9.

This product is subject to Antibody Design Labs Terms & Conditions of Sales available online at <http://www.abdesignlabs.com/terms/>.

© 2019 Antibody Design Labs. All rights reserved.